





Accessories for 2.813.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.813.0020

Qty.	Order no.	Description	
1 pcs	1.813.0010	813 Compact Autosampler Compact sample changer for 36 samples with a volume of 11 mL or 2.5 mL. Built-in peristaltic pump, 4 pre-programmed procedures, separate keyboard for making settings, control via «Remote» interface.	
1 pcs	6.1822.410	PEEK connection capillary Connection capillary for sample transfer with sample changers, 1 x M6, inner diameter 0.5 mm, for trace analysis	
2 pcs	6.1826.020	PVC pump tubing (blue/blue) For online IC instruments and automation in voltammetry.	
1 pcs	6.1831.050	PEEK capillary 0.5 mm i.d., 40 cm For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler	

1 pcs **6.1835.030** **PEEK sample tube**

For the Compact Autosampler.



1 pcs **6.2041.750** **Sample rack 36 x 11 mL**

Sample rack for 36 x 11 mL sample tubes, one-row



1 pcs **6.2621.060** **5/16 in. wrench**

For 5/16 in. screws. For IC instruments



2 pcs **6.2621.090** **1/2 in. wrench**

For 1/2 in. screws. For IC Sample Processors



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2743.057** **Sample tubes 11 mL**
200 pieces. For IC Sample Processors and VA Autosampler



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.030** **Coupling olive/UNF 10/32, 4x**
Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 pcs **6.2751.070** **Safety shield**

Safety shield for the 813 Compact IC Autosampler.



1 pcs

6.2752.010

Plug cover

Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor



0 pcs



6.2833.000

Needle/tube holder 1/8 in.

For IC/VA Sample Processors. For 1/8 in. tubes



Optional accessories

Order no.	Description	
2.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.	
2.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	
6.2141.170	Cable 797 VA Computrace to 813 Compact Sample Changer to 731 Relay Box	