






Accessories for 2.827.0114

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.827.0114

Qty.	Order no.	Description	
1 pcs	6.0228.010	Primatrode with NTC (fixed cable 1.2 m)	
		<p>Combined pH electrode with integrated temperature sensor (NTC) and fixed cable (1.2 m). This electrode is suitable for an introduction to GLP-compliant pH measurements in solutions that do not contain precipitates, proteins, or sulfides. This electrode is mechanically resistant thanks to the robust/unbreakable plastic shaft made of polypropylene and impact protection for the glass membrane. Reference electrolyte: c(KCl) = 3 mol/L, storage in storage solution.</p>	
1 pcs	6.2013.010	Clamping ring	
		For support rods with a diameter of 10 mm.	
1 pcs	6.2016.050	Support rod / 300 mm	
1 pcs	6.2021.020	Electrode holder	
		Electrode holder for 4 electrodes and 2 buret tips	
1 pcs	6.2161.010	Plug-in power supply unit EU 230V /6V DC	
		Power supply for 827 pH lab	


1 pcs

6.2621.070

5 mm hex key for IC Sample Processors



Optional accessories

Order no.	Description	
6.2307.230	Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit Maintenance kit for pH electrodes The kit contains: <ul style="list-style-type: none">• 50 mL cleaning solution• 50 mL 3M KCl solution• 50 mL storage solution• 2 Storage vessels• Instructions for use	