




Accessories for 2.833.0040

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.833.0040

Qty.	Order no.	Description	
7 pcs	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.	
1 pcs	6.1831.050	PEEK capillary 0.5 mm i.d., 40 cm For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler	
1 pcs	6.1831.060	PEEK capillary 0.5 mm i.d., 1 m For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	

1 pcs 6.2057.010 Holder for dialysis cell

For the dialysis cell in IC column chambers



1 pcs 6.2128.180 Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



1 pcs 6.2143.220 Cable 732/818 - Suppressor Module - Dialysis Unit

Connecting IC Detector - Suppressor Module and Dialysis Unit



1 pcs 6.2621.070 5 mm hex key for IC Sample Processors


1 pcs 6.2714.010 Dialysis membrane (cellulose acetate)

Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis



1 pcs 6.2729.100 Dialysis cell

With 4 x 6.2744.060 Blind stoppers. For inline dialysis



1 pcs 6.2744.000 Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



2 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



2 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.





4 pcs **6.2744.160** **Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



1 pcs **6.2831.010** **Plastic tweezers**
For dialysis and filtration membranes



Optional accessories

Order no.	Description	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2714.030	Dialysis membrane (polyamide) For inline dialysis and filtration	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	