





Accessories for 2.844.0120

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.844.0120

Qty.	Order no.	Description	
1 pcs	1.820.0500	820 IC Column Heating	
1 pcs	1.844.0120	844 UV/VIS Compact IC with post-column reactor and column heater Compact ion chromatograph for the photometric determination of anions and cations in the UV/VIS range. Very compact design. System with integrated dual-piston high-pressure pump, electrically operated injection valve and highly sensitive diode array detector. With completely metal-free flow path and additional column heating for two separation columns, PCR (Post-Column Reactor) and two-channel peristaltic pump.	
1 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32 For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	

1 pcs **6.1608.020** **Glass bottle / 1000 mL / GL 45**

Bottle for auxiliary solutions.



1 pcs **6.1608.070** **Eluent bottle / 2 L / GL 45**



1 pcs **6.1609.000** **Adsorbing tube, large and bent**

For filling with adsorber material.



1 pcs **6.1803.020** **PTFE capillary 0.97 mm i.d. / 5 m**

For all IC instruments



1 pcs **6.1803.050** **PTFE capillary 0.5 mm i.d. / 20 cm**

Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



1 pcs **6.1816.060** **Drain tubing for UV/VIS Compact IC and Bioscan**
Drainage tubing for draining off escaped liquid



2 pcs **6.1826.320** **Pump tubing LFL (orange/yellow) , 3-stopper**

For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



1 pcs **6.1831.010** **PEEK capillary 0.25 mm i.d. / 3 m**
For all IC components.



2 pcs **6.1831.100** **PEEK capillary 0.25 mm i.d., 1 m**



1 pcs **6.1831.110** **PEEK capillary 0.25 mm i.d., 30 cm**



1 pcs **6.1831.120** **PEEK capillary 0.25 mm i.d., 45 cm**



1 pcs **6.1834.010** **PTFE aspiration tube to IC Pump, 2.5 m**
For Compact IC, Personal IC, Basic IC, 6.5324.000 bottle holder.



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs **6.2027.030** **Column holder 8.5 mm**
For mounting columns in column chamber of IC instruments.



1 pcs **6.2027.040** **Column holder 11.3 mm**
For mounting columns in column chamber of IC instruments.



1 pcs **6.2057.050** **Holder for PCR-Module to 844 UV/VIS Compact IC**
Holding the Post Column Reactor in the UV/VIS Compact IC



1 pcs **6.2108.140** **Cable to column heating of 844**
Connecting cable 844 UV/VIS Compact IC to Column heater



1 pcs **6.2134.100** **RS-232 cable, 5 m, DB9 – DB9**
Connecting Metrohm RS-232 interface to Titrios, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



1 pcs **6.2617.010** **Tool for piston seal**
For removing and assembling the piston seal for all standard pump heads



1 pcs **6.2620.150** **Pulsation absorber (metal-free)**
For all ion chromatographs



1 pcs **6.2621.000** **Adjustable wrench**
Max. opening: 20 mm. For IC instruments



1 pcs **6.2621.030** **Hexagon key 4 mm**



2 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.030** **Coupling olive/UNF 10/32, 4x**
Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs 6.2744.180 Pump tubing connector with security lock and filter

For the connection of a pump tubing and a capillary with built-in filter



1 pcs 6.2816.020 Syringe 10 mL with Luer connection

For various applications in IC and VA



1 pcs 6.2821.090 Aspiration filter

Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



1 pcs 6.2821.130 Spare filter for inline filter

Spare filter for inline filter.



1 pcs 6.2836.000 Post column reactor

For mixing in coloring reagents in IC with UV/VIS detection.



1 pcs 6.2839.130 Flow cell 10 mm to UV/VIS instruments

PEEK flow cell

