





## Accessories for 2.844.0130

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.844.0130

Qty.	Order no.	Description	
1 pcs	1.820.0500	820 IC Column Heating	
1 pcs	1.844.0130	<b>844 UV/VIS Compact IC with column heater – SST</b> The 844 UV/VIS Compact IC with column heater - SST is a compact ion chromatograph for the photometric determination of polar substances, anions and cations in the UV/VIS range. The complete flow path is made of stainless steel and a column heater is installed. The 844 UV/VIS Compact IC with column heater - SST is equipped with a diode array detector.	
1 pcs	6.1602.160	<b>Bottle attachment GL 45 for eluents and reagents</b> For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
1 pcs	6.1608.070	<b>Eluent bottle / 2 L / GL 45</b>	

---

**1 pcs**      **6.1609.000**      **Adsorbing tube, large and bent**  
For filling with adsorber material.



---

**1 pcs**      **6.1816.060**      **Drain tubing for UV/VIS Compact IC and Bioscan**  
Drainage tubing for draining off escaped liquid



---

**1 pcs**      **6.1834.010**      **PTFE aspiration tube to IC Pump, 2.5 m**  
For Compact IC, Personal IC, Basic IC, 6.5324.000 bottle holder.



---

**1 pcs**      **6.2023.020**      **Clip for SGJ 14/15**  
Clip for SGJ 14/16



---

**1 pcs**      **6.2027.030**      **Column holder 8.5 mm**  
For mounting columns in column chamber of IC instruments.



---

**1 pcs**     **6.2027.040**     **Column holder 11.3 mm**  
For mounting columns in column chamber of IC instruments.



---

**1 pcs**     **6.2108.140**     **Cable to column heating of 844**  
Connecting cable 844 UV/VIS Compact IC to Column heater



---

**1 pcs**     **6.2134.100**     **RS-232 cable, 5 m, DB9 – DB9**  
Connecting Metrohm RS-232 interface to Titrios, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



---

**1 pcs**     **6.2617.010**     **Tool for piston seal**  
For removing and assembling the piston seal for all standard pump heads



---

**1 pcs**     **6.2620.000**     **Pressure screws (metal)**  
Pressure screws for ion chromatographs with stainless steel capillaries



---

**1 pcs**     **6.2620.160**     **Steel-Capillary 1/32 inch / 0.25 mm 30 cm**  
Capillary used with the stainless steel version of the UV/VIS Compact IC



---

**2 pcs**     **6.2620.170**     **Steel-Capillary 1/32 inch/0.25 mm 1 m**  
For capillary connections in the UV/VIS Compact IC



---

**1 pcs**     **6.2621.000**     **Adjustable wrench**  
Max. opening: 20 mm. For IC instruments



---

**1 pcs**     **6.2621.030**     **Hexagon key 4 mm**



---

**1 pcs**     **6.2621.050**     **1/4 in. wrench**  
For 1/4 in. screws. For IC instruments



---

**1 pcs**      **6.2744.270**      **Ferrule PEEK 1/32 inch 5 Pieces**  
For the connection of 1/32 inch steel capillaries



---

**1 pcs**      **6.2816.020**      **Syringe 10 mL with Luer connection**  
For various applications in IC and VA




---

**1 pcs**      **6.2821.090**      **Aspiration filter**  
Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



## Optional accessories

Order no.	Description	
<b>6.2621.040</b>	<b>Capillary Cutting Pincers</b>	
For steel capillaries. Used with IC instruments		