





Accessories for 2.881.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.881.0020

Qty.	Order no.	Description	
1 pcs	1.881.0020	881 Compact IC pro – Anion	 <p>The 881 Compact IC pro – Anion is an intelligent ion chromatograph of compact design with column oven for the determination of anions or polar substances with chemical suppression. It can be operated with various types of detection.</p>
2 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32	 <p>For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.</p>
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents	 <p>For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.</p>
1 pcs	6.1602.210	Attachment 5 x UNF 10/32 for «Waste Collector»	 <p>Attachment for holding up to 5 waste capillaries.</p>

2 pcs 6.1608.020 Glass bottle / 1000 mL / GL 45

Bottle for auxiliary solutions.



1 pcs 6.1608.070 Eluent bottle / 2 L / GL 45



1 pcs 6.1609.000 Adsorbing tube, large and bent

For filling with adsorber material.



1 pcs 6.1803.020 PTFE capillary 0.97 mm i.d. / 5 m

For all IC instruments



1 pcs 6.1803.040 PTFE capillary 0.5 mm i.d. / 1 m

Capillary for sample handling in IC.



1 pcs **6.1807.010** **Y-connector for tubing 6-9 mm i.d.**
Connector for waste tubing



1 pcs **6.1815.010** **Spiral band 0.5 m**
For holding together different cables or tubing.



2 pcs **6.1816.020** **Silicone tubing 6 mm i.d., 1 m**
For waste tubing.



1 pcs **6.1816.080** **Silicone tubing 6 mm i.d. / 2 m**

Drainage tubing to the "Waste Collector".



2 pcs **6.1826.420** **PharMed® pump tubing (orange/yellow), 3-stopper**

For MSM regeneration with peristaltic pump.



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs **6.2057.220** **Holder to Waste Collector**

Holder for fastening the Waste Collector to the IC instrument.



1 pcs **6.2151.020** **Cable USB A - USB B 1.8 m**
USB connecting cable



1 pcs **6.2251.000** **Colored sleeves for capillaries**
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 pcs **6.2617.010** **Tool for piston seal**
For removing and assembling the piston seal for all standard pump heads



2 pcs **6.2621.000** **Adjustable wrench**
Max. opening: 20 mm. For IC instruments



1 pcs **6.2621.030** **Hexagon key 4 mm**



1 pcs **6.2621.050** **1/4 in. wrench**
For 1/4 in. screws. For IC instruments



1 pcs **6.2621.080** **Capillary cutter**
For plastic capillaries. Used with IC instruments



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2626.000** **Drain nozzle (ProfIC)**
Drain nozzle for Professional IC instruments, for mounting on front panel



2 pcs **6.2739.000** **Wrench**
For tightening connectors



1 pcs **6.2743.080** **Stopper for overflow (ProfIC), 5 pieces**
For Professional IC instruments



1 pcs **6.2744.014** **Pressure screw 2x**
With UNF 10/32 connection. For the connection of PEEK capillaries



2 pcs **6.2744.020** **Luer/UNF coupling**
For IC instruments



1 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



2 pcs **6.2744.040** **2 x UNF 10/32 coupling**
For connecting 1/16 in. capillaries. For IC instruments



2 pcs **6.2744.070** **Pressure screw short**

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



2 pcs **6.2744.180** **Pump tubing connector with security lock and filter**
For the connection of a pump tubing and a capillary with built-in filter



1 pcs **6.2744.210** **Tubing adapter for aspiration filter (ProfIC)**
For Professional IC instruments



1 pcs **6.2765.000** **Container with GL 45 drainage to the "Waste Collector"**

Container with drainage to the "Waste Collector".



1 pcs **6.2816.020** **Syringe 10 mL with Luer connection**
For various applications in IC and VA



1 pcs **6.2816.040** **Purge tubing**
With PTFE tubing and Luer connection. For syringes. For aspirating eluents.









1 pcs **6.2821.090** **Aspiration filter**
Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



1 pcs **6.2821.130** **Spare filter for inline filter**
Spare filter for inline filter.



Optional accessories

Order no.	Description	
2.850.9010	IC Conductivity Detector	
	Compact and intelligent high performance conductivity detector for intelligent IC instruments. Outstanding temperature stability, the complete signal processing within the protected detector block and the latest generation of DSP – Digital Signal Processing – guarantee the highest precision of the measurement. No change of measuring ranges (not even automatic ones) is required, due to the dynamic working range.	
2.850.9110	IC Amperometric Detector	
	Compact and intelligent amperometric detector for intelligent IC instruments. Outstanding selectivity due to the four measuring modes: DC, PAD, flexIPAD and CV, as well as the excellent signal/noise ratio and the very fast stabilization of the measuring signal guarantee the highest in measurement precision.	
2.896.0010	896 Professional Detector – Conductivity	
	Intelligent stand-alone detector equipped with a high-performance IC Conductivity Detector. For use with intelligent IC instruments or as independent conductivity detector.	
2.896.0020	896 Professional Detector – Amperometry	
	Intelligent stand-alone detector equipped with the IC Amperometric Detector. With outstanding selectivity due to the four measuring modes: DC, PAD, flexIPAD and CV, as well as the excellent signal/noise ratio and very fast stabilization of the measuring signal guarantee highest precision of the measurement. For use with intelligent IC instruments or as independent detector.	
6.2617.040	Tool for piston seal, macro	
	For removing and assembling the piston seal for all macro pump heads	
6.2741.040	PE/PTFE piston seal, macro	
	For all macro pump heads	

6.2824.130 Macro pump head PEEK

Macro pump head for intelligent IC instruments, flow range 0.1...20 mL/min, maximum pressure 12.5 MPa.



6.5333.000 IQ/OQ Kit for IC

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.

