






Accessories for 2.917.2110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.917.2110

Qty.	Order no.	Description	
1 pcs	1.885.0010	885 Compact Oven Sample Changer The 885 Compact Oven Sample Changer is used for thermal sample preparation in Karl Fischer titration. In the oven the samples can be heated up to 250°C. The moisture contained in the sample evaporates and is transported by a dry carrier gas into the titration cell, where the analysis takes place.	
1 pcs	1.917.0010	917 Coulometer The 917 Coulometer is a compact titrator for coulometric water determination according to Karl Fischer. "Reduced to the max" – operating module, stirrer and pump are integrated in the 917 Coulometer for rapid reagent replacement. This saves a great deal of space on the laboratory table. Recommended measuring range: 10 µg...200 mg absolute water Coulometry is the ideal method for the determination of trace water in liquids, solids and gases. In addition, coulometry is an absolute method and thus no titer determination is necessary.	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0345.100	Generator electrode without diaphragm Generator electrode without diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for most applications where the sample is highly soluble in alcohol.	
1 pcs	6.1403.030	KF adsorber tube for coulometer cell For KF coulometers.	

1 pcs 6.1437.000 **SGJ stopper SGJ 14**



1 pcs 6.1439.010 **KF inlet and aspiration tube**
 KF inlet and aspiration tube for coulometer cell



2 pcs 6.1446.040 **Threaded stopper / M6**



1 pcs 6.1446.090 **Stopper / B-14/15 / M8**



2 pcs 6.1448.020 **Septum 16 mm, 5 pieces**
 Set of 5 pcs.



1 pcs 6.1448.067 **Screw cap with septum, 100 pcs.**
 Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces



1 pcs **6.1465.320** **KF titration vessel with 2 lateral openings / 80-250 mL / coulometric**
For 756 and 831 KF Coulometer



2 pcs **6.1602.105** **Bottle attachment / GL 45 / green**

For direct dosing from reagent bottles with GL 45 thread.



2 pcs **6.1602.145** **Drying flask cap**

Screw cap with GL 45 thread for the 6.1608.050 drying flask.



1 pcs **6.1608.023** **Amber glass bottle / 1000 mL / GL 45**
For Exchange Units. Bottle for auxiliary solutions.



1 pcs **6.1608.030** **Round glass bottle / 1000 mL / GL 45**



2 pcs **6.1608.050** **Drying bottle / 100 mL / GL 45**



1 pcs **6.1609.010** **Adsorber tube with tubing nipple**

Lid with tubing nipple.



1 pcs **6.1623.000** **Overflow protection for Ti Stand**

For 803 Ti Stand



1 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



1 pcs **6.1805.010** **FEP tubing / M6 / 13 cm**
With light and kink protection.



1 pcs **6.1805.060** **FEP tubing / M6 / 60 cm**
With light and kink protection



2 pcs **6.1805.080** **FEP tubing / M6 / 25 cm**
With light and kink protection.



2 pcs **6.1805.200** **PTFE tubing / M8 / 0.5 m**



1 pcs **6.1808.050** **Tubing adapter olive / M8 outer**
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 pcs **6.1819.030** **Tube**
For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0XX bottles



1 pcs **6.1819.050** **Tube**
For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



2 pcs **6.1821.040** **Filter tube**
Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer ovens.



1 pcs 6.1830.040 Heating tubing

Heating tubing with M6 thread for 885 Compact Oven Sample Processor in combination with a 915 Ti-Touch or 917 Coulometer.



1 pcs 6.1903.030 Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs 6.2023.020 Clip for SGJ 14/15

Clip for SGJ 14/16



2 pcs 6.2043.005 Holding clip for bottles

Holding spring for reagent bottles in exchange units.



1 pcs 6.2047.040 Titration vessel holder for 899/917

Titration vessel holder for coulometric cells for 899 Coulometer and 917 Ti-Touch



1 pcs 6.2049.050 Needle holder with Luer lock

Needle holder (73 mm) with Luer lock for the Karl Fischer ovens 860, 874 and 885.



1 pcs 6.2104.020 Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs 6.2104.120 Electrode cable / 1 m / H

With plug H. For connecting generator electrode - KF Coulometer.



1 pcs 6.2141.340 Remote cable 885 - Remote Box MSB

Remote cable for the connection of a 885 Compact Oven Sample Changer to Remote Box MSB (6.2148.010).



1 pcs 6.2148.010 Remote Box MSB

Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.



1 pcs 6.2164.010 Power supply unit 100 - 240 V/24 V DC

Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch



1 pcs 6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



1 pcs **6.2621.140** **Hexagon key 2.5 mm**



2 pcs **6.2701.040** **GL 18 screw cap with hole**

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



3 pcs **6.2713.000** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 14.



1 pcs **6.2713.010** **Sleeve**

Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.



1 pcs **6.2713.020** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



1 pcs **6.2724.010** **Dust filter**

Dust filter for Rancimats and Karl Fischer ovens.



1 pcs **6.2738.000** **Funnel**

For the KF coulometer.



1 pcs **6.2739.000** **Wrench**

For tightening connectors



1 pcs **6.2811.000** **Molecular sieve**

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



2 pcs **6.2816.030** **Needle with Luer connector**

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



1 pcs 6.2816.070 Septum piercing needle

Used with Karl Fischer ovens.



1 pcs 6.2816.080 Outlet needle

Used with KF Thermoprep and Oven Sample Processor.



2 pcs 6.2816.090 Syringe 5 mL

5 mL syringe made of PP, with Luer connector.



1 pcs 6.6066.000 USB flash drive for Ti-Touch

USB flash drive for 915 KF Ti-Touch and 916 Ti-Touch.



Optional accessories

Order no.	Description
-----------	-------------

2.142.0100 **Custom Q3X thermal printer**

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

6.2723.300 **Protective cover for 900 Touch Control**

Protective cover for 900 Touch Control

