




## Accessories for 6.00201.300

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.00201.300

Qty.	Order no.	Description	
1 pcs	6.1236.050	<b>Sleeve with SGJ 14/12 mm</b> Sleeve with SGJ 14/12 mm, polyethylene.	
1 pcs	6.2008.040	<b>Storage vessel</b>  Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.	

## Optional accessories

Order no.	Description	
6.02104.300	<b>Electrode cable plug-in head Q / plug P, 0.55 m</b>	
	Cable for connecting a dTrode (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	<b>Electrode cable plug-in head Q / plug P, 1.5 m</b>	
	Cable for connecting a dTrode (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	<b>Electrolyte 3 mol/L KCl (250 mL)</b>	
	Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)	
6.2323.000	<b>Storage solution</b>	
	Storage solution for all combined pH glass electrodes with reference electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$	
6.2325.000	<b>pHit kit</b>	
	Maintenance kit for pH electrodes	
	The kit contains:	
	<ul style="list-style-type: none"> <li>• 50 mL cleaning solution</li> <li>• 50 mL 3M KCl solution</li> <li>• 50 mL storage solution</li> <li>• 2 Storage vessels</li> <li>• Instructions for use</li> </ul>	

---

**6.2325.100****Cleaning solution 3 x 50 mL**

Reliable measuring results over long periods of time can be guaranteed only if the pH glass membrane and the diaphragm receive preventive and regular care. Cleaning by etching with toxic chemicals or applying mechanical treatment to the diaphragm is not only complicated and expensive, it also accelerates the ageing of the pH glass electrode.

The cleaning solution was developed for easy and gentle cleaning of pH glass electrodes. Regular use can considerably prolong their service life.

This cleaning solution is also part of the pHit kit.

