





Accessories for 6.00403.300

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.00403.300

Qty.	Order no.	Description	
1 pcs	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	
1 pcs	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.	

Optional accessories

Order no.	Description	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m	
	Cable for connecting a dTrode (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL)	
	Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	