







Accessories for 6.01400.100

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.01403.000	Titration head 6xNS14 / 3xNS9 (P&P)	
	Titration head with 6 NS14 and 3 NS 9 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place	
6.01900.010	Stirring propeller 30 mm ETFE	
	Stirring propeller 30 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	
6.02041.010	OMNIS sample rack 9 x 250 mL	
	OMNIS sample rack 9 x 250 mL for use with 250 mL glass sample beaker (P&P) or 250 mL PP sample beaker (P&P)	
6.02601.010	Gripper fingers 43 - 65 mm	
	Gripper fingers for OMNIS Sample Robot Pick&Place for gripping sample beakers with an outer diameter of 43 to 65 mm.	
6.02601.020	Gripper fingers 50 - 72 mm	
	Gripper fingers for OMNIS Sample Robot Pick&Place for gripping sample beakers with an outer diameter of 50 to 72 mm.	
6.02710.050	Dis-Cover lid for OMNIS 250 mL sample beaker, 9 pieces	
	Dis-Cover lid for 250 mL 6.01400.000, 6.01400.003 and 6.01400.100 sample beakers in the OMNIS Sample Robot Pick&Place system.	