







Accessories for 6.01400.300

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.01403.030	Titration head 3xNS14 / 4x6.4 mm (P&P)	
	Titration head with 3 NS14 and 4x 6.4 mm openings for using 120 mL beakers in Pick&Place modules of the OMNIS Sample Robot Pick&Place with two retracted dosing tips. Prepared for being equipped with rinsing and aspiration tubing.	
6.01404.030	Beaker adapter for OMNIS 120 mL sample beaker	
	Beaker adapter for centering the 120 mL sample beaker in the slide of the Pick&Place modules	
6.01900.030	Stirring propeller 20 mm ETFE	
	Stirring propeller, 20 mm, with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	
6.02041.030	OMNIS sample rack, 16 x 120 mL, (PP)	
	OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 16 sample beakers. The following sample beakers can be used: 6.01400.200, 6.01400.300, 6.01400.303. Plastic: Polypropylene (PP)	
6.02601.010	Gripper fingers 43 - 65 mm	
	Gripper fingers for OMNIS Sample Robot Pick&Place for gripping sample beakers with an outer diameter of 43 to 65 mm.	
6.02710.010	KF Dis-Cover lid for OMNIS 120 mL sample beaker, 16 pieces	
	Lid with seals for volumetric Karl Fischer Titration in 6.01400.200, 6.01400.300 and 6.01400.303 sample beakers in the OMNIS Sample Robot Pick&Place system.	

6.02710.040 Dis-Cover lid for OMNIS 120 mL sample beaker, 16 pieces

Dis-Cover lid for 120 mL 6.01400.200, 6.01400.300 and 6.01400.303 sample beakers in the OMNIS Sample Robot Pick&Place system.

