

## Accessories for 6.03001.120

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.03001.120

Qty.	Order no.	Description	
1 pcs	6.01503.120	<b>OMNIS 2 mL cylinder unit without accessories</b>	
		Intelligent 2 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.	
1 pcs	6.1446.030	<b>Ball stopper / B-14/(15)</b>	
		For the movable holding of buret tips in titration vessels	
2 pcs	6.1446.040	<b>Threaded stopper / M6</b>	
1 pcs	6.1543.200	<b>Titration tip / M6 thread</b>	
		Titration tip with antidiffusion valve.	

---

**2 pcs**     **6.1805.100**     **FEP tubing / M6 / 40 cm**  
With light and kink protection



---

**1 pcs**     **6.2008.030**     **Buret tip storage vessel for Dosing Unit**  
Together with 6.2008.050 buret tip holder provides a fastening for the buret tips on 6.3032.XXX Dosing units.



---

**1 pcs**     **6.2709.070**     **Guiding sleeve**  
For fixation of buret and dosing tips in titration heads and electrode holders



---

**1 pcs**     **6.2739.000**     **Wrench**  
For tightening connectors



---

**1 pcs**     **6.2803.010**     **Grease (2 g)**  
Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



## Optional accessories

Order no.	Description
6.05700.140	<b>Consumable Kit for OMNIS cylinder unit</b>

FEP M6 (40 cm) tubing connections and antidiffusion tip for OMNIS 6.03001.xxx cylinder units.

