






Accessories for 6.05800.030

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05800.030

Qty.	Order no.	Description	
3 pcs	6.01805.000	FEP tubing connection / M6 / 66.5 cm With light protection and protection against kinks.	
1 pcs	6.01805.010	PTFE tubing connection / M8 / 66.5 cm With light protection and protection against kinks.	
1 pcs	6.05001.010	"Peristaltic" (2-channel) pump module accessory set The accessory set contains all of the components required for the tube connections between a rinsing pump and an aspiration pump at the "Peristaltic" pump module. In addition to preassembled tubing with Luer connectors, each of these also includes two press clamps and pump tubing.	
1 pcs	6.1543.170	Aspiration tip / M8 thread Aspiration tip for the Sample Processors.	
2 pcs	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	

2 pcs **6.1621.100** **Cap for 10-L container**
Cap for 10-L container. For easy and safe transportation of filled containers.



2 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.



1 pcs **6.2709.070** **Guiding sleeve**
For fixation of buret and dosing tips in titration heads and electrode holders



3 pcs **6.2740.020** **Spray nozzle**
For the fine-spraying of the rinsing solution.

