





Accessories for 6.1204.180

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.2103.110	Contact pin M3, 2 mm	
	For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.	
6.2802.000	Polishing set for solid-state electrodes	
	Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)	
6.2802.020	Polishing set for Ultra Trace Electrode Tip	
	Polishing set to renew the surface of the Ultra Trace Electrode tip 6.1204.180	
6.2802.030	Polishing tape to Polishing Set for Ultra Trace Electrode	
	Polishing tape to Polishing Set for Ultra Trace Electrode. 12 µm particle size. used in the polishing cylinder. Set of 50	