

Accessories for 6.1245.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.0728.010	Ag/AgCl reference electrode without electrolyte	
Inner reference system without electrolyte for use with electrolyte vessel made of plastic 6.1245.010.		
6.0728.110	Ag/AgCl reference electrode without electrolyte	
Inner reference system without electrolyte for use with electrolyte vessel made of plastic 6.1245.010.		
6.0728.120	Ag/AgCl reference electrode	
Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.		
6.0728.130	LL-Ag/AgCl reference electrode for CVS	
Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.		