





Accessories for 6.1256.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1256.010

Qty.	Order no.	Description	
1 pcs	6.1244.050	Drive belts for rotating disk electrodes For Professional VA/CVS instruments (3 pcs.).	
1 pcs	6.1829.030	PTFE tubing / M6 / 21 cm, tapered	
1 pcs	6.2069.000	Park station for the measuring head of Professional VA/CVS instruments	
1 pcs	6.2709.110	Stopper for the electrode opening on Professional VA/CVS instruments Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).	

Optional accessories

Order no.	Description
6.1454.130 for KF and VA instruments	Form seal made of FPM for 61414xxx 