





## Accessories for 6.1462.070

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
6.1458.010	<b>Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings</b> With 6xSGJ 14 and 3x SGJ 9 apertures.	
6.1458.020	<b>Titration head /4 x M10 connectors</b> With 4 M10 screw nipples, including stoppers.	
6.1458.030	<b>KF titration head /4 x M10 connectors</b> With 4 M10 screw nipples, including stoppers.	
6.1458.040	<b>Titration Head, 3 x SGJ 14</b> Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.	
6.1458.070	<b>Titration head / 2 x SGJ 14</b> Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.	