



Accessories for 6.1462.240

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1808.220	Luer adapter for transfer head	
	Adapter with outside thread M6 and luer connection for 6.1462.240 Transfer head for using luer tools on multi-row racks	
6.1808.250	Pipette adapter 1-10 mL for transfer head	
	Adapter with outside thread M6 and Eppendorf conus to use with transfer head 6.1462.240 for usage of 10 mL Eppendorf tips	