



## Accessories for 6.1537.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.1537.010

Qty.	Order no.	Description	
1 pcs	6.1446.030	<b>Ball stopper / B-14/(15)</b> For the movable holding of buret tips in titration vessels	
1 pcs	6.1543.060	<b>Tip / M6 thread</b>  Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.	
1 pcs	6.1805.100	<b>FEP tubing / M6 / 40 cm</b> With light and kink protection	