



## Accessories for 6.1608.070

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.1602.150</b>	<b>Bottle attachment / GL 45 - 3 x UNF 10/32</b>	
	For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
<b>6.1602.160</b>	<b>Bottle attachment GL 45 for eluents and reagents</b>	
	For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	