


Accessories for 6.5104.040

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5104.040

Qty.	Order no.	Description	
1 pcs	6.2119.000	Mains adapter US 115V/9V AC With US mains plug. Input: 110...117 V AC. Output: 9 V AC/0.2 A. Used with 6.5104.040 differential amplifier with plug F / 115 V US	
1 pcs	6.2129.010	Differential amplifier, plug F Differential amplifier with plug F, for Titrandos. Power supply to be ordered seperately 6.2117.000 Power supply 230 V 6.2119.000 Power supply 115 V	