





Accessories for 6.5303.030

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5303.030

Qty.	Order no.	Description	
1 pcs	6.1803.000	Micro capillary tubing 1m To electrochemical detectors, 690 Ion Chromatograph, 698 Autosampler, 750 Autosampler	
1 pcs	6.2704.010	Seal for VA detector cell PTFE. Used with 656 Electrochemical Detector and 791 VA Detector	
1 pcs	6.2704.020	Seal for VA detector cell Used with 656 Electrochemical Detector and 791 VA Detector	

Optional accessories

Order no.	Description	
6.0727.000 ELCD reference electrode	Mini reference electrode for 791	
6.0805.010 Glassy carbon electrode	Mini Glassy Carbon Electrode for 791	
6.1204.130 Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	Silver Electrode Tip	
6.1204.140 Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	Gold Electrode Tip	