





## Accessories for 6.5610.020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5610.020

Qty.	Order no.	Description	
1 pcs	6.0340.000	<b>Double Pt-wire electrode for sample changer (fixed cable)</b>  Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.	
1 pcs	6.1458.030	<b>KF titration head /4 x M10 connectors</b>  With 4 M10 screw nipples, including stoppers.	
1 pcs	6.1805.060	<b>FEP tubing / M6 / 60 cm</b>  With light and kink protection	
1 pcs	6.1808.000	<b>Connection bushing M6</b>	

---

**1 pcs**     **6.1808.090**     **Thread adapter M8 outer / M6 inner**  
Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



---

**1 pcs**     **6.1821.000**     **Aspiration tube**  
For KF Sampler



---

**24 pcs**     **6.2037.040**     **Foil holder**  
Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer



---

**1 pcs**     **6.2820.000**     **Al foil**  
1000 circular sheets of 80 mm diameter. Thickness: 0.010 mm. Used together with 6.2037.030 Foil holder for 6.1432.210 Sample beaker. For 674 KF Sample Changer

