



## Accessories for 6.5622.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5622.000

Qty.	Order no.	Description	
1 pcs	6.1826.160	<b>Viton pump tubing</b>  For peristaltic pumps.	
1 pcs	6.2141.300	<b>Remote cable 843 Pump Station to Sample Processor</b> Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors	
1 pcs	6.2841.010	<b>External rinsing station for robotic arm with 6.1458.040 titration head</b>  The external rinsing station can be fitted on the side next to the Robotic Sample Processor tower, making it possible for external cleaning of the titration equipment to be carried out during additional actions on the sample rack.	