

## Accessories for 6.5702.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5702.010

Qty.	Order no.	Description	
1 pcs	6.1543.200	<b>Titration tip / M6 thread</b>  Titration tip with antidiffusion valve.	
2 pcs	6.1805.010	<b>FEP tubing / M6 / 13 cm</b>  With light and kink protection.	
2 pcs	6.1805.080	<b>FEP tubing / M6 / 25 cm</b>  With light and kink protection.	

---

**2 pcs**      **6.1805.100**      **FEP tubing / M6 / 40 cm**  
With light and kink protection



---

**1 pcs**      **6.1829.010**      **FEP aspiration tube M6, 0.25 m**  
With M 6 thread. For Dosing Unit



---

**4 pcs**      **6.2726.090**      **Antidiffusion valve**  
For 6.1541.XXX and 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200 antidiffusion buret tips

