

## Accessories for AUT.RRDE

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery AUT.RRDE

| Qty.  | Order no. | Description  |
|-------|-----------|--|
| 1 pcs | RRDE.PTPT | <b>Rotating Ring Disk Electrode, Platinum disk Platinum ring</b> |

The RRDE electrode tip consist of a 5 mm disc of platinum with a concentric platinum ring at a distance of 375 $\mu$ m, leading to a theoretical collection efficiency of 24.9%.



## Optional accessories

| Order no. | Description |
|-----------|-------------|
|-----------|-------------|

**RRDE.AUPT**      **Rotating ring disk electrode, Gold disk Platinum ring**

The RRDE electrode tip consist of a 5 mm disc of gold with a concentric platinum ring at a distance of 375 $\mu$ m, leading to a theoretical collection efficiency of 24.9%.



---

**RRDE.GCPT**      **Rotating ring disk electrode, Glassy Carbon disk Platinum ring**

The RRDE electrode tip consist of a 5 mm disc of Glassy Carbon with a concentric platinum ring at a distance of 375 $\mu$ m, leading to a theoretical collection efficiency of 24.9%.

