

Accessories for MVA-01


Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery MVA-01






Qty.	Order no.	Description
1 pcs	2.797.0010	797 VA Computrace for trace analysis

The 797 VA Computrace is a modern voltammetric measuring stand that is connected to a PC via a USB port. The computer software provided controls the measurement, records the measured data and evaluates it. Operation is most straightforward due to the well-laid-out structure of the program. The integrated potentiostat with galvanostat guarantees the highest sensitivity with reduced noise.

Voltammetry system for trace analysis and education. Complete accessories with VA Computrace software and all electrodes for a complete measurement system: Multi-Mode Electrode pro (MME pro), Ag/AgCl reference electrode and Pt auxiliary electrode.



Optional accessories

Order no.	Description	
2.800.0010	800 Dosino	
	Drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).	
2.843.0040	843 VA Membrane Pump Station for Compact VA systems	
	The 843 Membrane Pump Station has two built-in membrane pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace.	
	The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.	
2.843.0140	843 VA Peristaltic Pump Station for Compact VA systems	
	The 843 Peristaltic Pump Station has two built-in peristaltic pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace.	
	The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.	
2.846.0010	846 Dosing Interface	
	USB-capable control unit for connecting a maximum of four 800 Dosinos or 805 Dosimats for dosing and Liquid Handling tasks. A Touch Control or the connection to a PC with tiamo ™, MagIC Net, viva or 797 VA Computrace is required for operation.	
6.2055.100	Bottle holder, double	
	Bottle holder for 2 x 1 L bottles. For 6.5617.000 aspiration equipment	

6.3032.120 Dosing Unit 2 mL

Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.



6.5327.000 MVA-Hg, equipment for mercury determination

Complete set of accessories for the determination of mercury as per Application Bulletin 96. Contains rotating gold working electrode, reference electrode, glassy-carbon auxiliary electrode and other accessories.



6.5327.020 MVA-CVS, equipment for the CVS for 797 VA Computrace

Complete accessories set for the determination of organic additives in electroplating baths with CVS (Cyclic Voltammetric Stripping). Contains rotating platinum working electrode, reference electrode, platinum auxiliary electrode and additional accessories.

