






Accessories for FLWCL-PEEK

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
250BT	Screen-Printed Gold Electrode (Aux.:Pt; Ref.:Ag) / Ink BT	
	Screen-Printed Gold Electrode (Aux.:Pt; Ref.:Ag) / Ink BT. Suitable for working with microvolumes, for decentralized assays or to develop specific sensors.	
710	Screen-Printed Prussian Blue/Carbon Electrode	
	Screen-Printed Prussian Blue/Carbon Electrode. Ideal for the determination of hydrogen peroxide at a low detection potential.	
C223AT	Screen-Printed Gold Electrode (Aux.:Au; Ref.:Ag) d. 1.6 / Ink AT	
	Screen-Printed Gold Electrode (Aux.:Au; Ref.:Ag) / Ink AT. Working Electrode of 1.6 mm diameter.	
CAC	Cable connector for Screen-Printed Electrodes	
	Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.	
CAST	μStat Cable connector for Screen-Printed Electrodes	
	Connects any DropSens equipment with DropSens Screen-Printed Electrodes	