





Accessories for ProfIC Vario 3 Anion

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery ProfIC Vario 3 Anion

| Qty. | Order no. | Description | |
|-------|------------|---|--|
| 1 pcs | 2.850.9010 | IC Conductivity Detector |  <p>Compact and intelligent high performance conductivity detector for intelligent IC instruments. Outstanding temperature stability, the complete signal processing within the protected detector block and the latest generation of DSP – Digital Signal Processing – guarantee the highest precision of the measurement. No change of measuring ranges (not even automatic ones) is required, due to the dynamic working range.</p> |
| 1 pcs | 2.858.0020 | 858 Professional Sample Processor – Pump |  <p>The 858 Professional Sample Processor – Pump processes samples from 500 µL to 500 mL. The sample transfer takes place either with the installed bidirectional two-channel peristaltic pump or with an 800 Dosino.</p> |
| 1 pcs | 2.940.1510 | 940 Professional IC Vario ONE/SeS/PP/Prep 1 |  <p>The 940 Professional IC Vario ONE/SeS/PP/Prep 1 is the intelligent IC instrument with sequential suppression and peristaltic pump for suppressor regeneration in combination with Metrohm Inline Sample Preparation, e.g., Inline Ultrafiltration or Inline Dialysis. The instrument can be used with any separation and detection methods.</p> <p>Typical areas of application:</p> <ul style="list-style-type: none"> Anion or cation determinations with sequential suppression after Inline Ultrafiltration or Inline Dialysis |
| 1 pcs | 6.5330.100 | IC equipment: Inline dialysis |  <p>Accessory set for Inline Dialysis. For use with the 858 Professional Sample Processor and an additional 2-channel peristaltic pump.</p> |

Optional accessories

| Order no. | Description | |
|-------------------|--|---|
| 2.800.0010 | 800 Dosino |  |
| | Drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm). | |
| 2.941.0010 | 941 Eluent Production Module |  |
| | The 941 Eluent Production Module enables the automatic production of an eluent. It enables continuous working without manual intervention and guarantees stable retention times. It can be combined with all Metrohm IC instruments and is monitored by MagIC Net. | |
| 6.2041.440 | Sample rack 148 x 11 mL + 3 x 300 mL |  |
| | Sample rack for 148 x 11 mL samples with 3 rinsing beakers. | |
| 6.2041.760 | Sample rack 54 x 11 mL + 1 x 300 mL |  |
| | Sample rack with small diameter for 54 x 11 mL sample vessels in two rows and additionally 1 x 300 mL bottle. | |
| 6.2741.000 | Piston seal |  |
| | For all IC pumps. Particularly well-suited for eluents with organic solvents. | |
| 6.2743.050 | Sample tubes 11 mL |  |
| | 2000 pieces. For IC Sample Processors and VA Autosampler | |

6.2743.070 Stopper with perforation

2000 pieces. For IC Sample Processors



6.2744.080 M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



6.2832.000 MSM Rotor A

Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)



6.2842.000 MSM-HC Rotor A

Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)



6.2842.020 Adapter sleeve for Suppressor Vario

The use of the adapter sleeve makes it possible to install both the Standard MSM and the MSM-LC rotor in the MSM-HC housing.



6.2842.200**MSM-HC Rotor C**

Suppressor rotor for all IC instruments with MSM-HC drive (Metrohm Suppressor Module with high capacity) for cation suppression.



6.2844.000**MSM-LC Rotor A**

Suppressor rotor for use with 2 mm columns.



6.5330.130**IC equipment: Liquid Handling Station, left**

Accessory set for assembling the Liquid Handling Station on the Professional Sample Processor. The rinsing water supply can be accomplished by means of a peristaltic pump or a Dosino.

