





Accessories for 2.136.0730

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.136.0730

Qty.	Order no.	Description	
1 pcs	1.136.0700	Combustion Oven	
1 pcs	1.136.0720	LPG/GSS Module The LPG/GSS Module is a sample injection system for either liquefied gases or gases under enhanced pressure. The module is controlled automatically and moves the respective amounts of gas to the combustion oven for combustion. The Analytik Jena product has been optimized specially for coupling with ion chromatographs.	
1 pcs	6.7304.030	Oven coupling for LPG/GSS Connector for connecting the LPG/GSS module to the combustion oven. For Metrohm Combustion IC.	

Optional accessories

Order no.	Description	
6.7301.000	Combustion tube for Combustion IC	
Quartz combustion tube for the Combustion Module.		
