


Accessories for 6.5330.120

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5330.120

Qty.	Order no.	Description	
1 pcs	6.1562.130	Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.	
1 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32 For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
2 pcs	6.1608.070	Eluent bottle / 2 L / GL 45	
1 pcs	6.1803.070	PTFE capillary 0.97 mm i.d. / 40 cm Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit	

2 pcs **6.1803.080** **PTFE capillary 0.97 mm i.d. / 1 m**
Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



1 pcs **6.1819.110** **FEP aspiration tubing to canister**
Aspiration tubing to canister attachment with M6 connector. For the production of eluents.



1 pcs **6.1826.380** **Pump tubing LFL (gray/gray), 3 stoppers**
For Inline Dilution



1 pcs **6.1903.000** **Stirring bar / 8 mm**
Stirring Bar with magnetic core, PTFE coated, length 8 mm



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



1 pcs **6.2744.160** **Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



1 pcs **6.2744.200** **Adapter UNF 10/32 outer / M6 inner**
Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



1 pcs **6.2841.120** **Liquid Handling Station, left**
Combined rinsing and dilution station for sample changers with Swing Head. For assembly to the left of the tower.



1 pcs **6.3032.210** **Dosing Unit 10 mL**

Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.

