





Accessories for 2.756.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.756.0010

Qty.	Order no.	Description	
1 pcs	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.	
1 pcs	1.756.0010	756 KF Coulometer KF coulometer with built-in printer	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0344.100	Generator electrode with diaphragm Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for: <ul style="list-style-type: none"> • samples with very low absolute water contents < 50 ppm • coulometric water content determination in oil samples • applications where ketone reagents are used • applications where low-conducting reagents are used (addition of chloroform, xylene >10%) 	

1 pcs 6.1403.030 KF adsorber tube for coulometer cell

For KF coulometers.



1 pcs 6.1437.000 SGJ stopper SGJ 14



1 pcs 6.1446.060 Stopper / B-14/15 / M10








2 pcs 6.1448.020 Septum 16 mm, 5 pieces

Set of 5 pcs.



1 pcs 6.1464.320 KF titration vessel / 80-250 mL / coulometric



1 pcs	6.1903.030	Stirring bar / 25 mm	
		Stirring bar with magnetic core, PTFE covering, length 25 mm.	
1 pcs	6.2001.050	Support stand without support rod	
		Used with 756, 831 KF Coulometer and 728 Magnetic Swing-out Stirrer	
1 pcs	6.2013.010	Clamping ring	
		For support rods with a diameter of 10 mm.	
1 pcs	6.2016.030	Support rod / 250 mm	
1 pcs	6.2047.020	Titration vessel holder	
		Titration vessel holder for coulometric cells.Used with 756, 831 KF Coulometer	

1 pcs **6.2104.020** **Electrode cable / 1 m / F**
 For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs **6.2104.120** **Electrode cable / 1 m / H**
 With plug H. For connecting generator electrode - KF Coulometer.



1 pcs **6.2130.040** **Keyboard for 756**
 Keyboard for 756 KF Coulometer



3 pcs **6.2237.020** **Thermal paper 35 m**
 Sufficient for approx. 10'000 printed lines. For 756 Coulometer



1 pcs **6.2241.030** **Paper reel for thermal paper**



1 pcs **6.2701.040** **GL 18 screw cap with hole**
 Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



3 pcs **6.2713.000** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 14.



1 pcs **6.2713.010** **Sleeve**

Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.



1 pcs **6.2713.020** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



1 pcs **6.2730.030** **Stopper**

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 pcs **6.2738.000** **Funnel**

For the KF coulometer.



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



1 pcs 6.2816.030 Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



2 pcs 6.2816.090 Syringe 5 mL

5 mL syringe made of PP, with Luer connector.



Optional accessories

Order no.	Description
-----------	-------------

2.700.0020 **700 Dosino**

The 700 Dosino is a Dosing Unit drive. With fixed cable and mini-DIN plug.



2.703.0010 **703 Ti Stand with stirrer and pump**

Compact titration stand with magnetic stirrer and built-in pump for addition of solvents and for aspirating off titration vessel contents.



2.860.0010 **860 KF Thermoprep**

The 860 KF Thermoprep is designed for thermal sample preparation in Karl Fischer titration.

Numerous substances cannot be analysed by direct Karl Fischer titration as they are not soluble, react with the Karl Fischer reagent or release their water only very slowly or not until heated to high temperatures. Using tightly sealed sample vials, the sample is inserted in the oven. The samples can be subsequently analysed by any volumetric or coulometric KF titrator.



6.1439.010 **KF inlet and aspiration tube**

KF inlet and aspiration tube for coulometer cell



6.1464.323 KF titration vessel / 80-250 mL / coulometric

6.1465.320 KF titration vessel with 2 lateral openings / 80-250 mL / coulometric

For 756 and 831 KF Coulometer



6.2134.050 Connecting cable RS-232 DB9 to printer (DB25)

Connecting Metrohm RS-232 interface to e.g. Epson LX 300 printer. Used with Titrinos, and 756, 831 KF Coulometer



6.2134.100 RS-232 cable, 5 m, DB9 – DB9

Connecting Metrohm RS-232 interface to Titrinos, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



6.2148.000 Remote Box

For connection of a barcode reader and a PC keyboard to the KF Coulometer and to Titrinos



6.5406.000 Coulometric KF equipment, including generator electrode with diaphragm

Complete accessory set for KF coulometer.



6.5406.010 Coulometric KF equipment, including generator electrode without diaphragm

Complete accessory set for KF coulometer.



6.5617.000 Equipment for reagent replacement

Equipment for automatic reagent replacement with KF coulometer and connected Dosino.

