

Accessories for BWT-84000953

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery BWT-84000953

Qty.	Order no.	Description	
1 pcs	BWT-84000102	Laser safety goggles (1064 nm) Laser safety goggles for lasers with excitation at 1,064 nm – suitable for Class IIIB laser systems from B&W Tek.	
1 pcs	BWT-84000294	Battery Charging Cradle Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
1 pcs	BWT-84000295	Vials for liquids, 15 mm in diameter Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	BWT-84000296	Stylus/pen Stylus/pen combo for touch screens.	
1 pcs	BWT-84000298	NanoRam Carrying Case NanoRam hard plastic rugged carrying case.	
1 pcs	BWT-84000329	Pouch for 15 mm Vials Pouch for 15 mm vials.	
1 pcs	BWT-84000342	Right-Angle Adapter Right-angle adapter for hands-free measurement using handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	

1 pcs	BWT-84000343	Point-and-Shoot-Adapter	Point-and-shoot adapter for handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	
1 pcs	BWT-84000344	Vial Holder Adapter	Vial holder adapter, compatible with 15 mm-diameter vials.	
1 pcs	BWT-84000348	Polystyrene Validation Cap	Polystyrene Validation cap with integrated ASTM standard polystyrene.	
1 pcs	BWT-84000354	Pouch for adapters	Pouch for NanoRam adapters.	
1 pcs	BWT-84000403	10.8 V rechargeable lithium ion battery	Rechargeable lithium-ion battery, for use with NanoRam-1064, TacticID-GP Plus, TacticID-N Plus, TacticID-1064.	
1 pcs	BWT-84000414	18 V Power Adapter	Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.	
1 pcs	BWT-84000433	O-ring replacement kit for NanoRam instruments	NanoRam O-ring replacement kit – includes: O-rings, ID 9.00 mm, x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.	

1 pcs **BWT-84000858** **USB cable**

USB cable supporting micro USB port up to type A USB 2.0.



1 pcs **BWT-84000969** **Contact Immersion Probe (1064 nm)**

Handheld Raman Spectrometer 1064 accessories: Contact immersion probe, 6" length, working distance approx. 1 mm, operating pressure up to 2 bar at room temperature.



1 pcs **BWT-84000970** **NanoRam Large Bottle Adapter (1064)**

Handheld Raman Spectrometer 1064 accessories: Large bottle adapter with 10 mm working distance.





1 pcs **BWT-84000971** **NanoRam-1064 Software and Product Manual**

NanoRam-1064 accessories: USB flash drive containing NID EX software installation package; NanoRam-1064 Product Manual (PDF); NID EX Installation Instructions (PDF).



Optional accessories

Order no.	Description	
BWT-840000353	Shaft Protection Cap	
Handheld Raman shaft protective cap, package of 10.		
BWT-840000361	Handheld Laser Power Meter	
Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.		