





Accessories for BWT-840000972

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery BWT-840000972

Qty.	Order no.	Description	
1 pcs	BWT-840000101	Laser safety goggles (785 nm)	
		Laser safety goggles for 785 nm lasers – suitable only for Class IIIB laser systems from B&W Tek.	
1 pcs	BWT-840000295	Vials for liquids, 15 mm in diameter	
		Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	BWT-840000329	Pouch for 15 mm Vials	
		Pouch for 15 mm vials.	
1 pcs	BWT-840000403	10.8 V rechargeable lithium ion battery	
		Rechargeable lithium-ion battery, for use with NanoRam-1064, TacticID-GP Plus, TacticID-N Plus, TacticID-1064.	
1 pcs	BWT-840000406	TacticID Vial-Holder Adapter	
		TacticID Vial-Holder Adapter, compatible with 15 mm-diameter vials.	
1 pcs	BWT-840000414	18 V Power Adapter	
		Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.	

1 pcs	BWT- 840000416	TacticID Point-and-Shoot adapter	
		TacticID Point-and-Shoot adapter.	
1 pcs	BWT- 840000417	TacticID Validation Cap	
		TacticID Validation Cap with integrated ASTM-standard polystyrene. For use with TacticID-N Plus and TacticID-GP Plus.	
1 pcs	BWT- 840000418	TacticID-GP Protection Cap	
		TacticID-GP shaft protective cap.	
1 pcs	BWT- 840000422	TacticID Carrying Case	
		TacticID hard plastic rugged carrying case.	
1 pcs	BWT- 840000491	TacticID Right-Angle Adapter	
		TacticID Right-Angle Adapter.	
1 pcs	BWT- 840000858	USB cable	
		USB cable supporting micro USB port up to type A USB 2.0.	
1 pcs	BWT- 840000861	Micro USB flash drive	
		TacticID micro USB flash drive.	

1 pcs	BWT-840000912	O-ring replacement kit for TacticID instruments	
<p>TacticID O-ring replacement kit – includes: O-rings, ID 9.00 mm x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.</p>			
1 pcs	BWT-840000948	USB OTG Host Cable	
<p>USB OTG host cable for data transmission between a USB flash drive and a handheld system.</p>			
1 pcs	BWT-840000974	Complete Raman Spectral Library (TacticID Plus)	
<p>Complete Raman Spectral Library for TacticID Plus with > 10,000 spectra of common chemicals, narcotics, solvents and plastics. Includes NFPA704 and GHS safety instructions.</p>			
1 pcs	BWT-840000975	Common Chemicals, Narcotics, Explosives, Solvents and Plastics Spectral Library (TacticID Plus)	
<p>Raman Spectral Library for TacticID Plus with more than 1,200 spectra of common chemicals, narcotics, explosives, solvents and plastics. Includes NFPA704 and GHS safety instructions. May be subject to export control.</p>			
1 pcs	BWT-840000980	TacticID Plus Manuals	
<p>TacticID Plus Product Manual Set including Quick Reference Guide and TID Installation Instructions.</p>			
1 pcs	BWT-840000981	TacticID Plus Software	
<p>USB flash drive containing TacticID Plus TID software for PC, TID installation package and TacticID Plus Product Manuals.</p>			

Optional accessories

Order no.	Description	
BWT-840000294	Battery Charging Cradle Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
BWT-840000353	Shaft Protection Cap Handheld Raman shaft protective cap, package of 10.	
BWT-840000361	Handheld Laser Power Meter Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	
BWT-840000428	Wireless Router Wireless, battery-powered portable router; 3G/4G connectivity.	
BWT-840000498	Disposable Battery Cartridge TactiID battery cartridge for use with six type 123A disposable batteries (e.g. Surefire). Batteries not included.	
BWT-840000524	TacPac™ Adapter SERS analysis adapter for use with TacPac™-P SERS substrates.	
BWT-840000525	TacPac™ SERS substrate (package of 10) Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 10 units.	
BWT-840000526	TacPac™ SERS substrate (package of 30) Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 30 units.	

BWT-840000527 TacPac™ sample preparation kit (package of 10)

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 10.



BWT-840000528 TacPac™ Kit for TacticID

SERS analysis kit for use with Raman handheld spectrometers from B&W Tek. For Raman signal amplification and fluorescence reduction of hard-to-measure samples. For use with 785 nm laser excitation. Includes: tweezers, TacPac adapter, 10 SERS strips.



BWT-840000585 TacPac™ sample preparation kit (package of 30)

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 30.



BWT-840000603 Pouch for TacticID Adapters

TacticID (785 nm) Pouch for Adapters



BWT-840000620 TacPac™ SERS substrate (package of 100)

Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 100 units.



BWT-840000621 TacPac™ sample preparation kit (package of 100)

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 100.



BWT-840000669 TacticID Car Charger

TacticID mobile car charger.



BWT-840000938 MicroUSB-USB Flash Drive

Micro USB to USB flash drive for data storage.



BWT-840000976 Common Chemicals, Narcotics, Solvents and Plastics Spectral Library (TacticID Plus)

Raman Spectral Library for TacticID Plus with more than 1,180 spectra of common chemicals, narcotics, solvents and plastics. Includes NFPA704 and GHS safety instructions.



BWT-840000977 Raman Explosives Spectral Library (TacticID Plus)

Raman Explosives Spectral Library TacticID Plus including > 180 spectra. Includes NFPA704 and GHS safety instructions.

Note: This item may be subject to Export control to some territories.



BWT-840000978 Raman Gemstone Spectral Library (TacticID Plus)

Raman Gemstone Spectral Library for TacticID Plus including 467 spectra.



BWT-840000979 Raman Spectral Library of Chemical Warfare Agents (TacticID Plus)

Raman Spectral Library of Chemical Warfare Agents for TacticID Plus. Includes NFPA704 and GHS safety instructions. Note: This item may be subject to Export control to some territories.



BWT-840000982 TacticID Plus Remote Operation

Tablet with Android operating system pre-loaded with the TOS Plus Remote application installed for remote data acquisition and access. Includes camera and GPS functions. TID-WPR, Wireless router included.

