





## Accessories for 2.1001.2220

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.1001.2220

Qty.	Order no.	Description	
1 pcs	1.1001.0020	<b>OMNIS Titrator with stirrer</b>	
		OMNIS Titrator with magnetic stirrer, without function license	
1 pcs	6.00201.300	<b>dEcotrode plus</b>	
		Digital, combined pH electrode for OMNIS. This electrode is suitable for aqueous acid/base titrations. The fixed ground-joint diaphragm is insensitive to contamination. Reference electrolyte: c(KCl) = 3 mol/L, storage in storage solution.	
1 pcs	6.01600.010	<b>OMNIS Liquid Adapter</b>	
		Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	
1 pcs	6.01601.000	<b>Bottle cap, multi-use</b>	
		Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.	

---

**1 pcs    6.02004.010    Bottle holder**

Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.




---

**1 pcs    6.02005.000    OMNIS electrode holder**

Electrode holder for one OMNIS titration system.




---

**1 pcs    6.02005.010    Electrode Parking Station**

Electrode Parking Station for electrode storage on an OMNIS titration system.




---

**1 pcs    6.02100.010    Digital measuring module**

Digital measurement channel for one OMNIS Titrator or Titration Module for the connection of digital "dTrodes" electrodes.




---

**1 pcs    6.02104.300    Electrode cable plug-in head Q / plug P, 0.55 m**

Cable for connecting a dTrode (plug-in head Q) to the digital measuring module from OMNIS.




---

**1 pcs    6.02107.000    Network cable, 1 m**

Network cable (1 m) for connecting one OMNIS Sample Robot or OMNIS Titrator to an Ethernet Router, Ethernet Switch or local network



---

**1 pcs    6.03001.220    OMNIS 20 mL cylinder unit**

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.




---

**1 pcs    6.06000.000    OMNIS Software storage device**

License-based PC program for controlling complex OMNIS analysis systems.

- Sample-oriented operating procedures with numerous templates
- Graphical method editor with numerous templates
- Clearly organized and professional sample lists
- Management of instruments, work systems, sensors, solutions
- Dialog languages: German, English, Chinese

OMNIS  
A WHOLE NEW LEVEL OF PERFORMANCE

---

**1 pcs    6.06001.910    Function license Advanced**

"Advanced" function license for one OMNIS Titrator

- Liquid Handling (LQH)
- MEAS U/T/pH
- Calibration (Cal)
- Endpoint titration without conditioning (SET)
- Monotonic equivalence point titration (MET)
- Dynamic equivalence point titration (DET)

OMNIS  
A WHOLE NEW LEVEL OF PERFORMANCE

---

**1 pcs    6.06003.010    OMNIS Stand-Alone license**

Enables stand-alone operation of the OMNIS software on a Windows™ computer.

Features:

- The license already includes one OMNIS instrument license.
- Must be activated via the Metrohm licensing portal.
- Not transferable to another computer.

OMNIS  
A WHOLE NEW LEVEL OF PERFORMANCE

---

**1 pcs    6.1546.030    Extractor rod**

For PTFE piston of dosing unit



---

**1 pcs**     **6.1608.020**     **Glass bottle / 1000 mL / GL 45**

Bottle for auxiliary solutions.



---

**1 pcs**     **6.1903.020**     **Stirring bar / 16 mm**

Stirring Bar with magnetic core, PTFE coated, length 16 mm



---

**1 pcs**     **6.1903.030**     **Stirring bar / 25 mm**

Stirring bar with magnetic core, PTFE covering, length 25 mm.



---

**1 pcs**     **6.2013.010**     **Clamping ring**

For support rods with a diameter of 10 mm.



---

**1 pcs**      **6.2016.050**      **Support rod / 300 mm**



---

**1 pcs**      **6.2307.230**      **Buffer solutions pH 4, 7 and 9**

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



---

**1 pcs**      **6.2325.000**      **pHit kit**




Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use



## Optional accessories

Order no.	Description	
<b>6.05609.000</b>	<b>KFT equipment for OMNIS Titrator</b>	
	Complete accessory set for volumetric Karl Fischer titration.	
<b>6.05609.010</b>	<b>KFT equipment for OMNIS Titrator with Polytron</b>	
	Complete accessory set for volumetric Karl Fischer titration with Polytron.	
<b>6.06001.160</b>	<b>Function license KFT with conditioning</b>	
	Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.	