


















Accessories for 2.1001.3220

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1001.3220

Qty.	Order no.	Description	
1 pcs	1.1001.0020	OMNIS Titrator with stirrer	
		OMNIS Titrator with magnetic stirrer, without function license	
1 pcs	1.1006.0020	Rod Stirrer Titration	
		Rod stirrer for connection to one OMNIS Titrator (cable length 1.1 m).	
1 pcs	6.00203.300	dSolvotrode	
		Digital, combined pH electrode for OMNIS for all non-aqueous acid/base titrations. The glass membrane is optimized for poorly conducting solutions, and thanks to the flexible ground-joint diaphragm, the electrode is well suited for contaminated samples.	
		This electrode can be used with non-aqueous reference electrolytes (lithium chloride or tetraethylammonium bromide).	
		Storage in the respective reference electrolyte.	
1 pcs	6.01600.010	OMNIS Liquid Adapter	
		Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	

1 pcs	6.01601.000	Bottle cap, multi-use	<p>Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.</p> 
1 pcs	6.01900.010	Stirring propeller 30 mm ETFE	<p>Stirring propeller 30 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"</p> 
1 pcs	6.02004.010	Bottle holder	<p>Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.</p> 
1 pcs	6.02005.000	OMNIS electrode holder	<p>Electrode holder for one OMNIS titration system.</p> 
1 pcs	6.02005.010	Electrode Parking Station	<p>Electrode Parking Station for electrode storage on an OMNIS titration system.</p> 
1 pcs	6.02100.010	Digital measuring module	<p>Digital measurement channel for one OMNIS Titrator or Titration Module for the connection of digital "dTrodes" electrodes.</p> 

1 pcs	6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m	
Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.			
1 pcs	6.02107.000	Network cable, 1 m	
Network cable (1 m) for connecting one OMNIS Sample Robot or OMNIS Titrator to an Ethernet Router, Ethernet Switch or local network			
1 pcs	6.03001.220	OMNIS 20 mL cylinder unit	
Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.			
1 pcs	6.06000.000	OMNIS Software storage device	
License-based PC program for controlling complex OMNIS analysis systems.			
<ul style="list-style-type: none"> • Sample-oriented operating procedures with numerous templates • Graphical method editor with numerous templates • Clearly organized and professional sample lists • Management of instruments, work systems, sensors, solutions • Dialog languages: German, English, Chinese 			
1 pcs	6.06001.910	Function license Advanced	
"Advanced" function license for one OMNIS Titrator			
<ul style="list-style-type: none"> • Liquid Handling (LQH) • MEAS U/T/pH • Calibration (Cal) • Endpoint titration without conditioning (SET) • Monotonic equivalence point titration (MET) • Dynamic equivalence point titration (DET) 			
1 pcs	6.06003.010	OMNIS Stand-Alone license	
Enables stand-alone operation of the OMNIS software on a Windows™ computer.			
Features:			
<ul style="list-style-type: none"> • The license already includes one OMNIS instrument license. • Must be activated via the Metrohm licensing portal. • Not transferable to another computer. 			

1 pcs **6.1546.030** **Extractor rod**
For PTFE piston of dosing unit



1 pcs **6.1608.020** **Glass bottle / 1000 mL / GL 45**
Bottle for auxiliary solutions.



1 pcs **6.1903.020** **Stirring bar / 16 mm**
Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 pcs **6.1903.030** **Stirring bar / 25 mm**
Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2013.010** **Clamping ring**
For support rods with a diameter of 10 mm.






1 pcs **6.2016.050** **Support rod / 300 mm**



1 pcs **6.2320.000** **TEABr 0.4 mol/L in ethylene glycol (250 mL)**
Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$



Optional accessories

Order no.	Description	
6.05609.000	KFT equipment for OMNIS Titrator	
	Complete accessory set for volumetric Karl Fischer titration.	
6.05609.010	KFT equipment for OMNIS Titrator with Polytron	
	Complete accessory set for volumetric Karl Fischer titration with Polytron.	
6.06001.160	Function license KFT with conditioning	
	Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.	