

## Accessories for 2.848.3010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.848.3010

Qty.	Order no.	Description
1 pcs	1.802.0010	<b>802 Rod Stirrer for 727 Ti Stand</b>

Rod stirrer with separate 6.2137.000 Speed controller and 6.1909.010 Propeller stirrer.



1 pcs	1.804.0010	<b>804 Ti Stand without stand rod</b>
-------	------------	---------------------------------------

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.



1 pcs	1.848.0010	<b>848 Titrino plus</b>
-------	------------	-------------------------

Compact titrator for

- dynamic (DET) titrations with automatic equivalence point finding
- monotonic (MET) titrations with automatic equivalence point finding
- endpoint titrations (SET)



---

**1 pcs    6.0229.010    Solvotrode easyClean (fixed cable 1.2 m)**

Combined pH electrode with flexible easyClean diaphragm and fixed cable (1.2 m) for all non-aqueous acid/base titrations. The glass membrane is optimized for poorly conducting solutions, and thanks to the easy-to-clean easyClean diaphragm, the electrode is well suited for heavily contaminated samples, such as used oil. This electrode can be used with non-aqueous reference electrolytes (lithium chloride or tetraethylammonium bromide). Storage in the respective reference electrolyte.




---

**1 pcs    6.1909.010    Stirring propeller / 96 mm**


Propeller stirrer, length from lower edge of SGJ: 96 mm. For use with beakers with 722, 802 Rod Stirrer




---

**1 pcs    6.2001.060    Support stand without support rod**

To attach an 801, 804, 803 Stirrer to a Titrande, Titrimo plus, Dosimat plus or 856 and 867.




---

**1 pcs    6.2013.010    Clamping ring**

For support rods with a diameter of 10 mm.




---

**1 pcs    6.2016.070    Support rod / 400 mm**


---

**1 pcs**     **6.2021.020**     **Electrode holder**  
Electrode holder for 4 electrodes and 2 buret tips



---

**1 pcs**     **6.2147.000**     **Numerical USB keypad**  
Numerical USB keypad for Titrino plus, Dosimat plus, 862 Compact Titrosampler, 869 Compact Sample Changer, 863 Compact Autosampler and 899 Coulometer. With two integrated USB connectors.



---

**1 pcs**     **6.2151.100**     **Adapter USB MINI (OTG) - USB A**  
For connecting USB instruments.



---

**1 pcs**     **6.2320.000**     **TEABr 0.4 mol/L in ethylene glycol (250 mL)**  
Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol),  $c(\text{TEABr}) = 0.4 \text{ mol/L}$



---

**1 pcs**     **6.2621.070**     **5 mm hex key for IC Sample Processors**



---

**1 pcs**      **6.2621.130**      **Hexagon key 2 mm**  
2 mm.



---

**1 pcs**      **6.2723.310**      **Protective cover for Titrino plus**

Protective cover for display and keypad for all instruments of the Titrino plus family.

Suitable for: 848 Titrino plus, 862 Compact Titrosampler, 863 Compact Autosampler, 865 Dosimat plus, 869 Compact Sample Changer, 870 KF Titrino plus, 876 Dosimat plus, 877 Titrino plus, 899 Coulometer.



---

**1 pcs**      **6.2739.010**      **Wrench**  
For exchange units.



---

**1 pcs**      **6.3026.220**      **Exchange Unit 20 mL**

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



---

**1 pcs**





**6.6061.010**

**Titrimo plus methods USB Stick**

Contains the FOOD PAC, WINE PAC and OIL PAC methods for direct import into the Titrimo plus.



## Optional accessories

Order no.	Description	
<b>2.142.0100</b>	<b>Custom Q3X thermal printer</b>	
	Compact printer with USB interface for	
	<ul style="list-style-type: none"> <li>• 900 Touch Control</li> <li>• 915 KF Ti-Touch</li> <li>• 916 Ti-Touch</li> <li>• 917 Coulometer</li> <li>• 877 / 848 Titrino plus</li> <li>• 865 / 876 Dosimat plus</li> <li>• 91X Meter (cable 6.2151.140)</li> <li>• Eco Dosimat / Titrator</li> <li>• 862 Compact Titrosampler</li> <li>• 870 KF Titrino plus</li> <li>• 899 Coulometer</li> </ul>	
	Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.	
<b>6.1115.000</b>	<b>Optrode</b>	
	Optical sensor for photometric titrations offering 8 different wavelengths. The wavelength can be switched using the software (tiamo 2.5 or higher) or with a magnet. The glass shaft is completely solvent-resistant and easy to clean. For example, this space-saving sensor is suitable for:	
	<ul style="list-style-type: none"> <li>• Non-aqueous titrations in accordance with USP or EP</li> <li>• Determinations of carboxyl end groups</li> <li>• TAN/TBN in accordance with ASTM D974</li> <li>• Sulfate determination</li> <li>• Fe, Al, Ca in cement</li> <li>• Water hardness</li> <li>• Chondroitin sulfate in accordance with USP</li> </ul>	
	The sensor is not suitable for determinations of concentrations via measurement of color intensity (colorimetry).	
<b>6.2104.020</b>	<b>Electrode cable / 1 m / F</b>	
	For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
<b>6.2141.230</b>	<b>Remote cable Compact Sample Changer/Compact Autosampler - Titrino plus/Pump Station</b>	
	Cable for the remote connection between Compact Sample Changer/Compact Autosampler and Titrino plus and/or Pump Station.	

---

**6.2141.240 Remote cable Dosimat plus - Compact Sample Changer/Titrino plus**

Cable for the remote connection between Dosimat plus and Compact Sample Changer and/or Titrino plus using 6.2141.260 Cable.



---

**6.2141.250 Remote cable Compact Sample Changer - Titrino**

Cable for the remote connection between Compact Sample Changer and Titrinos with 25-pin remote socket.



---

**6.2141.260 Remote cable Compact Sample Changer - Dosimat plus and Titrino plus**

Cable for remote connection between Compact Sample Changer and Dosimat plus and further devices (e.g. Titrino plus or Pump Station).



---

**6.2151.120 Cable USB MINI (OTG) - USB B, 1.8 m**

For connecting USB instruments.



---

**6.3026.110 Exchange Unit 1 mL**

Exchange Unit with integrated data chip with 1 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.



---

**6.3026.150 Exchange Unit 5 mL**

Exchange Unit with integrated data chip with 5 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



---

**6.3026.210 Exchange Unit 10 mL**

Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



---

**6.3026.250 Exchange unit 50 mL**

Exchange unit with integrated data chip with 50 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion tip and standard amber glass reagent bottle.

