






Accessories for 2.919.0120

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.919.0120

Qty.	Order no.	Description	
1 pcs	1.919.0020	919 IC Autosampler plus	 <p>The 919 IC Autosampler plus fulfills the requirements of laboratories with medium sample numbers. It enables automation of the full range of Metrohm IC instruments.</p>
1 pcs	6.1458.080	Transfer head for 919 IC Autosampler plus for VA and CVS applications.	
1 pcs	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m For all IC instruments	
1 pcs	6.1815.010	Spiral band 0.5 m For holding together different cables or tubing.	
2 pcs	6.1826.020	PVC pump tubing (blue/blue) For online IC instruments and automation in voltammetry.	

1 pcs **6.1831.020** **PEEK capillary 0.75 mm i.d., 22 cm**
For all IC components.



1 pcs **6.1835.050** **PEEK sample tube**
For sample changers for transferring larger-sized sample volumes, e.g. for brightener determination.



1 pcs **6.2041.490** **Sample rack VA 28 x 50 mL**
Sample rack for 28 samples in 50 mL sample vessels (6.2747.010), single-row.



1 pcs **6.2151.000** **Cable USB A – mini-DIN 8-pin**
Controller cable



1 pcs **6.2621.080** **Capillary cutter**
For plastic capillaries. Used with IC instruments



1 pcs **6.2621.090** **1/2 in. wrench**
For 1/2 in. screws. For IC Sample Processors



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2621.130** **Hexagon key 2 mm**
2 mm.



1 pcs **6.2739.000** **Wrench**
For tightening connectors



2 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



2 pcs 6.2744.060 Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



1 pcs 6.2744.080 M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



2 pcs 6.2744.160 Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 pcs 6.2744.290 T connector 3 × UNF 10/32 with holder



4 pcs 6.2747.010 Sample vessel 50 mL

With blue screw cap and mL graduation. Set of 25 pieces.



1 pcs 6.2752.110 Holder for T connector to 919 IC Autosampler plus



1 pcs **6.2833.030** **Needle holder 1/8 in. for 863 Compact Autosampler**
Needle holder 1/8 in. for 863 Compact Autosampler

