

Accessories for 2.863.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.863.0010

Qty.	Order no.	Description
------	-----------	-------------

1 pcs	1.863.0010	863 Compact IC Autosampler
-------	------------	----------------------------

The 863 Compact IC Autosampler is the ideal help for routine analysis. It enables automation of the full range of Metrohm IC instruments.



2 pcs	6.1826.360	Pump tubing LFL (white/white), 3 stoppers
-------	------------	---

For Sample Processors



1 pcs	6.1831.050	PEEK capillary 0.5 mm i.d., 40 cm
-------	------------	-----------------------------------

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler



1 pcs **6.1831.060** **PEEK capillary 0.5 mm i.d., 1 m**
For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



1 pcs **6.2151.100** **Adapter USB MINI (OTG) - USB A**
For connecting USB instruments.



1 pcs **6.2151.110** **Cable USB MINI (OTG) - USB A, 1.8 m**
For connecting USB instruments.



1 pcs **6.2621.090** **1/2 in. wrench**
For 1/2 in. screws. For IC Sample Processors



1 pcs **6.2621.130** **Hexagon key 2 mm**
2 mm.



1 pcs **6.2621.140** **Hexagon key 2.5 mm**



1 pcs **6.2739.000** **Wrench**
For tightening connectors



1 pcs **6.2743.057** **Sample tubes 11 mL**
200 pieces. For IC Sample Processors and VA Autosampler



1 pcs **6.2743.077** **Stopper with perforation**
200 pieces. For IC Sample Processors



1 pcs 6.2744.034 Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



1 pcs 6.2744.070 Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs 6.2744.160 Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 pcs 6.2833.030 Needle holder 1/8 in. for 863 Compact Autosampler

Needle holder 1/8 in. for 863 Compact Autosampler



1 pcs 6.2846.010 Sample needle made of zirconium oxide with PEEK tip

Sample needle for Sample Processors and Autosampler in ion chromatography. The PEEK tip guarantees secure piercing of the vial covering.



Optional accessories

Order no.	Description
-----------	-------------

6.1826.310	Pump tubing LFL (orange/green), 3-stopper
-------------------	--

Pump tubing for suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



6.1826.320	Pump tubing LFL (orange/yellow), 3-stopper
-------------------	---

For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



6.1826.330	Pump tubing LFL (orange/white), 3 stoppers
-------------------	---

For all IC instruments with peristaltic pump.



6.1826.340	Pump tubing LFL (black/black), 3 stoppers
-------------------	--

For the sample solution in inline dialysis.



6.1826.380 Pump tubing LFL (gray/gray), 3 stoppers

For Inline Dilution



6.1826.390 Pump tubing LFL (yellow/yellow), 3 stoppers

For the sample solution in Inline Ultrafiltration.



6.2141.230 Remote cable Compact Sample Changer/Compact Autosampler - Titrimo plus/Pump Station

Cable for the remote connection between Compact Sample Changer/Compact Autosampler and Titrimo plus and/or Pump Station.



6.2141.270 Cable Compact IC Autosampler - Remote Connection on IC Instruments

Enables remote control of the Compact IC Autosampler from IC Instruments



6.2151.020 Cable USB A - USB B 1.8 m

USB connecting cable



6.2629.000 Screw nut to needle holder IC

On Sample Processors



6.2723.310 Protective cover for Titrino plus

Protective cover for display and keypad for all instruments of the Titrino plus family.

Suitable for: 848 Titrino plus, 862 Compact Titr sampler, 863 Compact Autosampler, 865 Dosimat plus, 869 Compact Sample Changer, 870 KF Titrino plus, 876 Dosimat plus, 877 Titrino plus, 899 Coulometer.



6.2743.040 Sample tubes 2.5 mL

2.5 mL. 2000 pieces. For IC Sample Processors



6.2743.047 Sample tubes 2.5 mL

2.5 mL. 200 pieces. For IC Sample Processors



6.2743.050 Sample tubes 11 mL

2000 pieces. For IC Sample Processors and VA Autosampler

