

Accessories for 2.135.0205



Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.135.0205

Qty.	Order no.	Description
1 pcs	1.135.0205	Balance for 864 Robotic Balance Sample Processor, 115V The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 115V.



Optional accessories

Order no.	Description	
6.2067.100	Scale pan for sample rack 6.2068.010 The scale pan is designed for the 120 mL PP beaker (6.1459.300) and is placed into the sample rack of the 864 Balance Sample Processor.	
6.2068.010	Sample rack 20 x 120 mL for 864 Balance Sample Processor Sample rack for 864 Balance Sample Processor with 20 x 120 mL PP beakers (6.1459.300) for weighing and analysing of the samples, 20 pipetting tips (6.1562.240), 20 x scale pans (6.2067.100) and 20 additional 120 mL PP beakers for samples.	
6.2711.080	Drip pan for 864 Robotic Balance Sample Processor The drip pan guards Robotic Balance Sample Processor and balance against contamination / corrosion by chemicals.	