





Accessories for 2.776.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.776.0010

Qty.	Order no.	Description	
1 pcs	1.776.0010	776 Dosimat Universal dispensing unit for liquid handling in laboratories; also usable for titrators with pulse control, with automatic stopcock switching, including push-button cable for manual dispensing control, keypad and 6.3026.220 Exchange Unit 20 mL.	
1 pcs	6.2107.000	Push-button cable for Dosimat For manual operation of Dosimats	
1 pcs	6.2149.000	Keyboard for 765/776 Keyboard for 765, 776 Dosimat	
1 pcs	6.2739.010	Wrench For exchange units.	

1 pcs

6.3026.220

Exchange Unit 20 mL

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

