



Accessories for 6.0908.110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0908.110

Qty.	Order no.	Description	
1 pcs	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	
1 pcs	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.	

Optional accessories

Order no.	Description
-----------	-------------

6.2301.060 **KCl ion standard 250 mL**

KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25 °C).



6.2324.010 **Conductivity standard 100 µS/cm 250 mL**

Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.



6.2324.110 **Conductivity standard 100 µS/cm, 5 x 30 mL**

Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.

