


Accessories for 2.136.0800

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.136.0800

Qty.	Order no.	Description
1 pcs	1.136.0800	MMS 5000 Autosampler
		MMS 5000 Autosampler from Analytik Jena for use with the Metrohm Combustion IC for fully automatic analysis of liquid and solid samples. In order to match the modular Multi-Matrix sampler with the correct sample type, either the liquid kit (6.7303.000) or the solid kit (6.7302.000) must be used.



Optional accessories

Order no.	Description	
6.7301.150	Substrate made of quartz for liquid samples Substrate for liquid sample injections, 3 pcs.	
6.7302.000	Kit for solid samples for MMS 5000, 35 positions Kit for solid samples for MMS 5000 with rack for sample boats, grippers and 10 quartz glass boats.	
6.7302.100	Quartz boat 40 x 9 mm Quartz glass boat for sample injection in Combustion Module, 10 pcs.	
6.7302.310	Sample rack 35 positions for solid samples For solids.	
6.7303.000	Kit for liquid samples for MMS 5000, 112 positions Kit for liquid samples for MMS 5000 with rack for sample vessels, 100 µL syringe with replaceable needle and 100 snap-ring bottles.	
6.7303.100	Snap-ring bottle ND11, 1.5 mL For liquid holder with MMS 5000, 100 pcs.	
6.7303.110	Septum cap ND11 For liquid holder with MMS 5000, 100 pcs.	

6.7303.200 **50 μ L syringe with replaceable needle**

Replacement syringe for MMS 5000 liquid setup



6.7303.210 **100 μ L syringe with replaceable needle**

Replacement syringe for MMS 5000 liquid setup



6.7303.220 **Needles for 50 μ L/100 μ L syringes, 3 pcs.**

Replacement needle for syringe for MMS 5000 liquid setup



6.7303.310 **Sample rack 112 positions for liquid samples**

For liquid samples.

