

## Accessories for 2.885.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.885.0010

Qty.	Order no.	Description	
1 pcs	1.885.0010	<b>885 Compact Oven Sample Changer</b> The 885 Compact Oven Sample Changer is used for thermal sample preparation in Karl Fischer titration. In the oven the samples can be heated up to 250°C. The moisture contained in the sample evaporates and is transported by a dry carrier gas into the titration cell, where the analysis takes place.	
1 pcs	6.1446.170	<b>Heating tube stopper</b> for coulometric KF vessels in connection with heatable tubing connection	
1 pcs	6.1448.067	<b>Screw cap with septum, 100 pcs.</b> Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces	
2 pcs	6.1602.145	<b>Drying flask cap</b> Screw cap with GL 45 thread for the 6.1608.050 drying flask.	
2 pcs	6.1608.050	<b>Drying bottle / 100 mL / GL 45</b>	

---

**1 pcs**     **6.1805.010**     **FEP tubing / M6 / 13 cm**  
With light and kink protection.



---

**1 pcs**     **6.1805.060**     **FEP tubing / M6 / 60 cm**  
With light and kink protection



---

**2 pcs**     **6.1805.080**     **FEP tubing / M6 / 25 cm**  
With light and kink protection.



---

**1 pcs**     **6.1808.040**     **Thread adapter M6 outer / M8 inner**  
Outer thread M6, inner thread M8.



---

**1 pcs**     **6.1808.050**     **Tubing adapter olive / M8 outer**  
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



---

**2 pcs**     **6.1821.040**     **Filter tube**  
Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer ovens.



---

**1 pcs**     **6.1830.030**     **Heating tubing**  
Heating tubing with M6 thread for 860 KF Thermoprep, 885 Compact Oven Sample Processor and 874 USB Oven Sample Processor.



---

**2 pcs**     **6.2043.005**     **Holding clip for bottles**  
Holding spring for reagent bottles in exchange units.



---

**1 pcs**     **6.2049.050**     **Needle holder with Luer lock**  
Needle holder (73 mm) with Luer lock for the Karl Fischer ovens 860, 874 and 885.



---

**1 pcs**     **6.2151.100**     **Adapter USB MINI (OTG) - USB A**  
For connecting USB instruments.



---

**1 pcs**     **6.2420.007**     **Sample glass with thread, 8 mL, 100 pieces**  
Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



---

**1 pcs**     **6.2621.140**     **Hexagon key 2.5 mm**



---

**1 pcs**     **6.2724.010**     **Dust filter**  
Dust filter for Rancimats and Karl Fischer ovens.



---

**1 pcs**     **6.2739.000**     **Wrench**  
For tightening connectors



---

**1 pcs    6.2811.000    Molecular sieve**

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



---

**1 pcs    6.2816.070    Needle for gas inlet**

Used with KF Thermoprep.



---

**1 pcs    6.2816.080    Outlet needle**

Used with KF Thermoprep and Oven Sample Processor.

