





Accessories for 2.889.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.889.0020

Qty.	Order no.	Description	
1 pcs	1.889.0020	889 IC Sample Center – cool	
		The 889 IC Sample Center – cool is the appropriate automation solution when you have only a very small amount of sample. In comparison with the 889 IC Sample Center, it possesses in addition a cooling function and is thus the ideal sample changer for biochemically relevant or thermally unstable samples.	
1 pcs	6.1608.110	Wash bottle 250 mL to IC Sample Center	
		Bottle for washing solution in the IC Sample Center.	
2 pcs	6.1807.010	Y-connector for tubing 6-9 mm i.d.	
		Connector for waste tubing	
2 pcs	6.1816.070	Silicon tubing 8 mm i.d. /1 m	
		For the Sample Center	

1 pcs **6.1831.010** **PEEK capillary 0.25 mm i.d. / 3 m**
For all IC components.



2 pcs **6.2041.200** **Sample rack 48 x 11 mm to IC Sample Center**
Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.



1 pcs **6.2151.020** **Cable USB A - USB B 1.8 m**
USB connecting cable



1 pcs **6.2743.107** **PP Vials 0.7 mL to IC Sample Center, 100 pieces**
100 sample vials.



1 pcs **6.2743.127** **Septum to 6.2743.1XX Si/PTFE**
100 pieces



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2816.200** **Air needle 80 mm, black**
To IC Sample Center



Optional accessories

Order no.	Description	
6.1562.170	Tefzel transfer tubing 1000 μ L To IC Sample Center.	
6.1829.060	Capillary set for syringe unit For the IC Sample Center	
6.2743.100	PP Vials 0.7 mL to IC Sample Center, 1000 pieces 1000 sample vials.	
6.2743.110	PP Vials 0.3 mL to IC Sample Center, 1000 pieces 1000 sample vials.	
6.2743.117	PP Vials 0.3 mL to IC Sample Center, 100 pieces 100 sample vials.	
6.2743.120	Septum to 6.2743.1XX Si/PTFE 1000 pieces	

6.2816.210 Needle with capillary, 15 μ L

To IC Sample Center



6.5333.000 IQ/OQ Kit for IC

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.



6.7303.100 Snap-ring bottle ND11, 1.5 mL

For liquid holder with MMS 5000, 100 pcs.



6.7303.110 Septum cap ND11

For liquid holder with MMS 5000, 100 pcs.

