




Accessories for 2.942.0070

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.942.0070

Qty.	Order no.	Description	
1 pcs	1.942.0070	942 Extension Module Vario LQH	
		The 942 Extension Module Vario LQH is an extension module for the Professional IC Vario line of instruments. It expands the possibilities of Metrohm Inline Sample Preparation "MISP" and Liquid Handling and opens a Professional IC system for Online Monitoring.	
		Typical applications:	
		<ul style="list-style-type: none"> • Extended Inline Sample Preparation • Online monitoring • Dilution without sample changers • Fully automatic sample derivation reactions 	
1 pcs	6.1562.160	Transfer tubing 15 mL, 2 x M6	
		For liquid handling in ion chromatography. To avoid contamination of the Dosino.	
5 pcs	6.1608.100	PE bottle 50 mL	
		Bottle for Extension Module Liquid Handling.	
1 pcs	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m	
		For all IC instruments	
1 pcs	6.1825.230	PEEK sample loop 10 µL	
		For injection valve, with 2 PEEK pressure screws	

2 pcs 6.1826.340 Pump tubing LFL (black/black), 3 stoppers

For the sample solution in inline dialysis.



1 pcs 6.1831.010 PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



1 pcs 6.2057.200 Bottle holder to Extension Module - Liquid Handling



1 pcs 6.2156.060 Cable Extension Module - Professional IC, 40 cm

Cable connecting the extension module with a Professional IC instrument



1 pcs 6.2251.000 Colored sleeves for capillaries

Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 pcs 6.2744.010 Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs 6.2744.034 Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



2 pcs 6.2744.070 Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



2 pcs 6.2744.080 M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



6 pcs 6.2744.090 Pressure screw long

Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



2 pcs 6.2744.160 Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 pcs **6.2744.200** **Adapter UNF 10/32 outer / M6 inner**
Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)





1 pcs **6.2762.000** **Mixing vessel 50 mL to Extension Module - Liquid Handling**



1 pcs **6.2762.010** **Cover to mixing vessel 6.2762.000**



Optional accessories

Order no.	Description	
6.2103.170	Adapter to Leak Sensor at Professional IC Adapter to connect two leak sensors to a Professional IC instrument	
6.2156.070	Cable Extension Module - Professional IC, 1 m Cable connecting the extension module with the Professional IC instrument	
6.5332.000	IC tubing kit for second extension module	