


## Accessories for 6.2824.130

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.2824.130

Qty.	Order no.	Description	
2 pcs	6.2741.040	<b>PE/PTFE piston seal, macro</b> For all macro pump heads	
1 pcs	6.2824.160	<b>Outlet valve PEEK</b> Metal-free check valve for the outlet of IC pump heads	
1 pcs	6.2824.170	<b>Inlet valve PEEK</b> Metal-free check valve for the inlet of IC pump heads	