




Accessories for 2.843.0130

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.843.0130

Qty.	Order no.	Description	
1 pcs	1.843.0120	843 Pump Station (peristaltic) The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.	
1 pcs	6.1543.120	Aspiration tip / M8 thread Aspiration tip if 6.1805.200 siphoning tubing is used	
2 pcs	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	
2 pcs	6.1621.100	Cap for 10-L container Cap for 10-L container. For easy and safe transportation of filled containers.	

4 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



1 pcs **6.1805.200** **PTFE tubing / M8 / 0.5 m**



3 pcs **6.1805.420** **FEP tubing / M6 / 48 cm**

With light and kink protection



1 pcs **6.1808.240** **Distributor Rinsing / Aspiration for 869 Compact Sample Changer**

Distributor with 3 x M6, 1 x M8 and 2 x 1/8 connections for the use of rinsing and / or aspiration equipment on the 869 Compact Sample Changer.



2 pcs **6.1812.000** **PTFE tubing 4 mm / 6 mm / 4 m**



4 pcs **6.1820.050** **M6 / Olive adapter**

For 6.1826.100 pump tubing, to peristaltic pumps.



1 pcs **6.1826.100** **PP pump tubing 6.4 mm / 1 m**

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



2 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.






3 pcs **6.2740.020** **Spray nozzle**

For the fine-spraying of the rinsing solution.



Optional accessories

Order no.	Description	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus Cable for remote connection between Compact Sample Changer and Dosimat plus and further devices (e.g. Titrino plus or Pump Station).	
6.2141.300	Remote cable 843 Pump Station to Sample Processor Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors	
6.2740.030	Rinsing nozzle The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.	
6.2761.000	Splash protection for 843 Pump Station (Peristaltic) The peristaltic pump runs with mounted splash protection only.	