

Accessories for 2.843.0150

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.843.0150

Qty.	Order no.	Description	
1 pcs	1.843.0120	843 Pump Station (peristaltic) The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.	
1 pcs	6.1543.170	Aspiration tip / M8 thread Aspiration tip for the Sample Processors.	
2 pcs	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	
2 pcs	6.1621.100	Cap for 10-L container Cap for 10-L container. For easy and safe transportation of filled containers.	

4 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



3 pcs **6.1805.060** **FEP tubing / M6 / 60 cm**
With light and kink protection



3 pcs **6.1805.420** **FEP tubing / M6 / 48 cm**
With light and kink protection



1 pcs **6.1805.510** **PTFE tubing / M8 / 60 cm**

With protection against kinks.



1 pcs **6.1808.170** **Distributor for Sample Processors for rinsing and aspirating**

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



2 pcs **6.1812.000** **PTFE tubing 4 mm / 6 mm / 4 m**



4 pcs **6.1820.050** **M6 / Olive adapter**

For 6.1826.100 pump tubing, to peristaltic pumps.



1 pcs **6.1826.100** **PP pump tubing 6.4 mm / 1 m**

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



2 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.



1 pcs **6.2141.300** **Remote cable 843 Pump Station to Sample Processor**

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors




3 pcs **6.2740.020** **Spray nozzle**

For the fine-spraying of the rinsing solution.



Optional accessories

Order no.	Description	
6.2740.030	Rinsing nozzle	
<p>The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.</p>		
