






Accessories for 2.851.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.851.0020

Qty.	Order no.	Description	
1 pcs	1.801.0010	801 Stirrer	
		Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).	
1 pcs	1.851.0010	851 Titrando	
		Coulometry is the ideal method for water determination in liquids, solids and gases when it comes to water determination in the trace range (10 µg to 10 mg absolute water). In addition, coulometry is an absolute method and thus no titer determination is necessary.	
		Coulometric titrations are carried out easily and quickly with the 851 Titrando .	
		Recommended measuring range: 10 µg - 200 mg absolute water	
1 pcs	1.900.0010	900 Touch Control	
		Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick.	
		Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian.	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry	
		Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0344.100	Generator electrode with diaphragm	
		Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for:	
		<ul style="list-style-type: none"> • samples with very low absolute water contents < 50 ppm • coulometric water content determination in oil samples • applications where ketone reagents are used • applications where low-conducting reagents are used (addition of chloroform, xylene >10%) 	

1 pcs **6.1403.030** **KF adsorber tube for coulometer cell**

For KF coulometers.



1 pcs **6.1437.000** **SGJ stopper SGJ 14**



1 pcs **6.1446.060** **Stopper / B-14/15 / M10**



2 pcs **6.1448.020** **Septum 16 mm, 5 pieces**

Set of 5 pcs.



1 pcs **6.1464.320** **KF titration vessel / 80-250 mL / coulometric**




1 pcs **6.1903.030** **Stirring bar / 25 mm**

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2001.060** **Support stand without support rod**

To attach an 801, 804, 803 Stirrer to a Titrande, Titrino plus, Dosimat plus or 856 and 867.



1 pcs **6.2013.010** **Clamping ring**

For support rods with a diameter of 10 mm.



1 pcs **6.2016.030** **Support rod / 250 mm**



2 pcs **6.2043.005** **Holding clip for bottles**

Holding spring for reagent bottles in exchange units.



1 pcs **6.2047.020** **Titration vessel holder**
Titration vessel holder for coulometric cells. Used with 756, 831 KF Coulometer



1 pcs **6.2103.130** **Adapter red 2 mm plug / 4 mm socket**
For connecting plug B (4 mm) to 2 mm socket.



1 pcs **6.2103.140** **Adapter black 2 mm plug / B socket 4 mm**
For connecting plug B (4 mm) to 2 mm socket.



1 pcs **6.2104.020** **Electrode cable / 1 m / F**
For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs **6.2104.120** **Electrode cable / 1 m / H**
With plug H. For connecting generator electrode - KF Coulometer.



1 pcs **6.2701.040** **GL 18 screw cap with hole**
Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



3 pcs **6.2713.000** **Sleeve**
For the fat-free ground-joint connection with the standard ground-joint SGJ 14.



1 pcs **6.2713.010** **Sleeve**
Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.



1 pcs **6.2713.020** **Sleeve**
For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



1 pcs **6.2730.030** **Stopper**
Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 pcs **6.2738.000** **Funnel**
For the KF coulometer.



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



2 pcs 6.2816.030 Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



2 pcs 6.2816.090 Syringe 5 mL

5 mL syringe made of PP, with Luer connector.



1 pcs 6.6064.010 USB Memory Stick for 900 Touch Control